

LISS


Laparoscopic Irrigation Suction Set

Fairmont Medical's Laparoscopic Irrigation Suction Set is designed to meet all your minimally invasive surgical needs.

For over 25 years Fairmont Medical has been pioneering the innovative use of suction irrigation devices. Our extensive range of **LISS** sets are designed to provide cost effective, reliable, leak free laparoscopic surgical procedures.

LISS is known around the world for it's quality and consistent performance during minimal invasive procedures, such as laparoscopic cholecystectomy.

- Ergonomic trumpet valve accommodates even the most delicate hand sizes, ensuring comfortable use, minimising fatigue.
- Easy to feel, concave/convex colour coded buttons
- Unique detachable suction probe connection system offers fast change over probe choice
- Multi-use suction tube connection allows the use of other suction devices
- Atraumatic tip design with distal vents prevents wall occlusion
- Soft drapeable tubing
- Double spike bag available
- Single use, sterile, ready to use
- Phthalate (DEHP) and latex free
- Quality Australian designed and manufactured



LISS


"Respected around the world and exported to more than in 30 countries."





Fairmont
MEDICAL

LISS Laparoscopic Irrigation Suction Sets & Accessories


Permanently Connected Probe - Sterile, Disposable

Product	5mm Probe OD Product Code	10mm Probe OD Product Code	Probe Length	Bag Spike	Tubing Length	Box Qty
	LISS4021	-	380mm	Single	3.7m	10
	LISS4022	LISS4402	380mm	Double	3.7m	10
	LISS4052	-	500mm	Double	3.7m	10

Detached Probe - Sterile, Disposable

	LISS4121	-	380mm	Single	3.7m	10
	LISS4122	LISS4422	380mm	Double	3.7m	10
	LISS4127	-	440mm	Double	3.7m	10


No Probe - Sterile, Disposable

Product	Product Code	Bag Spike	Tubing Length	Box Qty
	LISS4011	Single	3.7m	10
	LISS4012	Double	3.7m	10

Irrigation Suction Probe

Product	Sterile, Disposable Product Code	Non-Sterile, Reusable Product Code	Length (without connector)	Probe OD	Box Qty
	-	LISS4007	300mm	5mm	1 each
	-	LISS4008	325mm	10mm	1 each
	-	LISS4002	400mm	5mm	1 each
	LISS4004	-	340mm	3mm	5
	LISS4009	-	340mm	5mm	20
	LISS4017	-	430mm	5mm	10
	LISS4117	-	340mm	10mm	5
	LISS4050	-	500mm	5mm	10
	LISS4417	-	430mm	10mm	10

SpeciJar - Collection Jar - Sterile, Disposable

Product	Product Code	Soft Tubing Length	Tube ID / OD	Fluid Container Volume	Box Qty
	DSJ5007	150mm	8mm/ 12.5mm	70ml	30

Fairmont
MEDICAL

Corporate Office
Fairmont Medical
11 Scoresby Road Bayswater 3153 VIC
Ph 1300 972 088
Ph International +613 9720 8840
Fx International +613 9237 9399

UK Head Office
Fairmont Medical Products Ltd
Elizabeth House
28 Baddow Road Chelmsford
Essex CM2 0DG
Ph: +44 (0)1245 206800

