



aircair O₂ Nasal Cannula


Fairmont Medical's Nasal Cannula delivers oxygen when a low flow, low / medium concentration is required.

Fairmont Medical's **Aircair** O₂ Nasal Cannula range is designed to provide a comfortable fit on all patients and are ideal for any oxygen therapy requiring up to 5-6 L/min flow rate.

- Designed for a comfortable fit on adults, children and infants
- Available in soft or standard tips with 2.1 metre oxygen tubing
- Soft nasal prongs designed for extra patient comfort
- Single use, ready to use
- Australian designed quality



Non-Sterile, Disposable

Product	Product Code	Description	Box Qty
	DNC9510	Disposable Adult Nasal Oxygen Cannula - Standard Tip - 2.1m O ₂ Tubing	100
	DNC9530	Disposable Adult Nasal Oxygen Cannula - Soft Tip - 2.1m O ₂ Tubing	100
	DNC9540	Disposable Paediatric Nasal Oxygen Cannula - Standard Tip - 2.1m O ₂ Tubing	100
	DNC9560	Disposable Paediatric Nasal Oxygen Cannula - Soft Tip - 2.1m O ₂ Tubing	100
	DNC9570	Disposable Infant Nasal Oxygen Cannula - Soft Tip - 2.1m O ₂ Tubing	100

Fairmont
MEDICAL

Corporate Office
Fairmont Medical
11 Scoresby Road Bayswater 3153 VIC
Ph 1300 972 088
Ph International +613 9720 8840
Fx International +613 9237 9399

UK Head Office
Fairmont Medical Products Ltd
Elizabeth House
28 Baddow Road Chelmsford
Essex CM2 0DG
Ph: +44 (0)1245 206800

