


Bile Bag

Fairmont Medical offers a quality cost effective Bile Bag to suit all your needs.

Fairmont Medical's Bile Bag provides a cost effective alternative to post operative biliary drainage collection. An included easy to manipulate T-tube adaptor assists in the emptying of bag contents.

- Transparent see through front panel
- Clear graduated markings
- Fitted with a multi-purpose tapered connector
- 500ml capacity
- In-built easy to attach holding strap eyelets
- Latex free
- Single use, sterile, ready to use
- Australian designed quality



Sterile, Disposable			
Product	Product Code	Capacity	Box Qty
	DBB7500	500mL	50



Corporate Office
Fairmont Medical
11 Scoresby Road Bayswater 3153 VIC
Ph 1300 972 088
Ph International +613 9720 8840
Fx International +613 9237 9399

UK Head Office
Fairmont Medical Products Ltd
Elizabeth House
28 Baddow Road Chelmsford
Essex CM2 0DG
Ph: +44 (0)1245 206800

