

bladavac

Disposable Bladder Evacuator

Fairmont Medical's Bladder Evacuator provides a simple solution for bladder flushing and tissue collection.

Fairmont Medical's **Bladavac** offers urological surgical teams a practical cost effective solution for the removal of tissue debris during TURP's, with full control. **Bladavac** comes complete with a filter trap system which prevents tissue debris returning to the surgical site. The soft easy to squeeze bottle gives the surgeon full control over this process.

- Easy to fit connection to resectoscope instrument
- Soft easy to squeeze bottle design allows fast and efficient removal of tissue and unwanted body fluids
- Easy to fit screw top cap
- Single use, sterile, ready to use
- Quality Australian designed

bladavac

"An effective way to flush the bladder and collect tissue samples."

"Prevents cross contamination and associated risks of exposing theatre personnel."



Sterile, Disposable		
Product	Product Code	Box Qty
	DBE5100	15

Fairmont
MEDICAL

Corporate Office
Fairmont Medical
11 Scoresby Road Bayswater 3153 VIC
Ph 1300 972 088
Ph International +613 9720 8840
Fx International +613 9237 9399

UK Head Office
Fairmont Medical Products Ltd
Elizabeth House
28 Baddow Road Chelmsford
Essex CM2 0DG
Ph: +44 (0)1245 206800

