

ABSORBEST



Fairmont
MEDICAL

ABSORBEST DryMax

No more slippery floors. No more soaked, cold and dripping cloths and sheets. May we introduce DryMax Operation, a range of super-absorbent mats and pads that efficiently manages excess fluids that occur on floors, beds and tables in operating theatres, emergency situations or during medical examinations.

DryMax super-absorbent products quickly absorb and encapsulate large amounts of water-based liquids, such as saline, urine and blood, including contaminated parts, without leaking. The result is a drier, safer and more hygienic working environment – with improved infection control – which also means less time and money spent on cleaning.


- Superior absorption – superior retention
- Comfortable, user-friendly – easy to apply
- No dripping – easy and safe to change and discard
- Fewer changes – less work and lower material costs
- Simplicity – easy to use

DryMax super-absorbent pads improve hygiene and infection control within operating rooms.




DryMax Sterile

Super-absorbent pads for use with the patient. Available in two sizes, with or without tape.




DryMax XL Medium

Anti-slip super absorbent floor mat for collection of splashing liquids. Contains a nonwoven top layer and holds up to 3L liquid.



DryMax XL Basic


Anti-slip floor mat for collection of splashing liquids.




DryMax 2.4

Super-absorbent floor mat with high retention capacity, in a handy size.

DryMax - Sterile, Disposable

Product	Product Code	Water Absorption	Saline Absorption Rate	Dimensions	Box Qty
	AF50243 Maxi with Tape	6L	2.4L	720mmx 370mm	80
	AF50163 Midi with Tape	4L	1.6L	560mmx 370mm	100

DryMax XL Range - Non-Sterile, Disposable

Product	Product Code	Water Absorption	Saline Absorption Rate	Dimensions	Box Qty
	AF22201 (Blue) DryMax XL Basic	1.5 L	1.5L of Fluid	1000mm x 750mm	50
	AF22203-1 DryMax XL Medium	3L	3L	1000mm x 750mm	20

DryMax 2.4 - Non Sterile, Disposable

Product	Product Code	Water Absorption	Saline Absorption Rate	Dimensions	Box Qty
	AF20241 (Green)	9L	9L of Water 3.5L of Saline	720mm x 370mm	100

