

Capnography Monitoring


Fairmont Medical's Capnography Masks provide breath-to-breath ventilation data.



Capnography Mask Range

Fairmont Medical's Capnography Masks range provides breath to breath monitoring of nose and mouth breathing. The CO₂ mask is a safe consistent means of monitoring breathing of a non-intubated patient, allowing sampling of exhaled carbon dioxide from point of exhalation, while at the same time administering oxygen.











- Extensive range of mask configurations available
- All encompassing nose and mouth monitoring approach
- Extension tubing allows for easy ISO luer lock connection
- See through design allows clear patient visualisation
- Soft conforming shape allows for patient comfort
- Patient strap locking, prevents pull outs
- Filter option available
- Phthalate (DEHP) and latex free
- Quality Australian designed and manufactured



Adult Cut Down Mask



Capnography Monitoring

Super Cut Down Masks - Non-Sterile, Disposable				Full Face Masks - Non-Sterile, Disposable			
Product	Product Code	Description	Box Qty	Product	Product Code	Description	Box Qty
	DAM6180	Disposable Adult Super Cut Down Capnography mask - 3m O ₂ tubing & CO ₂ sample line outlet	60		DAM9080	Disposable Adult Full Face Capnography Mask - 3m O ₂ tubing & CO ₂ sample line outlet	30
	DAM6182	Disposable Adult Super Cut Down Capnography Mask - 3m CO ₂ sample line & 3m O ₂ tubing	60		DAM9081	Disposable Paediatric Full Face Capnography Mask - CO ₂ sample line outlet & 3m O ₂ tubing	30
	DAM6183	Disposable Adult Super Cut Down Capnography Mask - CO ₂ sample line outlet, filter & 3m O ₂ tubing	60		DAM9082	Disposable Adult Full Face Capnography Mask - 3m CO ₂ sample line & 3m O ₂ tubing	30
	DAM6186	Disposable Adult Super Cut Down Capnography Mask - 3m CO ₂ sample line, filter & 3m O ₂ tubing	60		DAM9083	Disposable Adult Full Face Capnography Mask - CO ₂ sample line outlet, filter & 3m O ₂ tubing	30
					DAM9086	Disposable Adult Full Face Capnography Mask - 3m CO ₂ sample line, filter & 3m O ₂ tubing	30
					DAM9096	Disposable Paediatric Full Face Capnography Mask - 3m CO ₂ sample line, filter & 3m O ₂ tubing	30

