

Capnography Monitoring

Fairmont Medical's Capnography Nasal Cannula provides flexibility for patient CO₂ monitoring.



Capnography Nasal Cannula

Introducing Fairmont Medical's Nasal Cannula which provides CO₂ monitoring and oxygen delivery allowing greater patient access.

- Patient comfort fit via soft nasal prongs
- Compact dual patient delivery and monitoring
- Unique oxygen wall divide giving accurate CO₂ monitoring
- Easy ISO luer lock connector
- Flexible choice of configurations on offer
- Filter option available
- Phthalate (DEHP) and latex free



CO₂ Sampling Lines


Fairmont Medical's extensive range of CO₂ Sampling Lines are available in both PVC and PE lined tubing configurations. This provides choice for those seeking cost effective capnography monitoring.



- Maximum accuracy when it counts via PE lined tubing
- Compatible with popular monitoring models utilising our ISO luer lock connectors
- Longer lengths enables monitoring during MRI use
- Filter option available
- Phthalate (DEHP) and latex free



Capnography Monitoring

Capnography Nasal Cannula - Non-Sterile, Disposable					
Product	Product Code Male Luer Lock	Product Code Female Luer Lock	Description	Length	Box Qty
	DNC9600	DNC9610	Disposable Adult with CO ₂ Monitoring Nasal Oxygen Cannula	3.1 metres	80
	DNC9601	DNC9611	Disposable Adult with CO ₂ Monitoring Nasal Oxygen Cannula and Filter	3.1 metres	80

CO ₂ Sample Lines - Non-Sterile, Disposable					
Product	Product Code	Tubing	Length	Connectors	Box Qty
	DAF9006	Single Lumen PVC	3m	2 Male Luer Locks	100
	DAF9007	Dual Lumen, Outer Wall PVC, Inner Wall PE	3m	2 Male Luer Locks	100
	DAF9009	Dual Lumen, Outer Wall PVC, Inner Wall PE	9m	2 Male Luer Locks	50
	DAF9011	Single Lumen PVC	3m	1 Male Luer Lock 1 Female Lock	100


CO ₂ Sample Lines with Filter - Non-Sterile, Disposable					
Product	Product Code	Tubing	Length	Connectors	Box Qty
	DAF9008	Dual Lumen, Outer Wall PVC Inner Wall PE	150mm	1 Male Luer Lock 1 Female Luer Lock	100
	DAF9106	Single Lumen	3m	2 Male Luer Locks	100

Mixing Cannula

Mixing Cannula is easily connected to standard luer slip syringes. The semi-rigid tube is used to aspirate pharmaceutical fluids from vials and ampoules during the preparation of medicines.

- Cost effective
- Supplied sterile ready to use
- Universal luer slip connection
- Quality Australian designed and manufactured



Sterile, Disposable			
Product	Product Code	Length	Box Qty
	MIX1001	140mm	400

