

dasflo Disposable Arthroscopy Sets

For over 25 years Fairmont Medical has been providing orthopaedic surgeons with a safe cost effective method of maintaining a clear view during major joint surgery without the need for an expensive mechanical pump.

- Soft easy to squeeze pump chamber system
- Allows for quick high flow irrigation delivery and rapid refill through it's uniquely designed one-way valve system
- Easy instrument connection and flexible use made possible by silicone tube
- Optional double bag spike allows high volume irrigation bag set up
- Single use, sterile, ready to use
- Phthalate (DEHP) and latex free
- Quality Australian designed and manufactured

Fairmont Medical's Disposable Arthroscopy Sets provide a cost effective solution and are efficient to use eliminating the need for a mechanical pump.








dasflo *"Saves valuable surgical operating room time."*





Disposable Arthroscopy Sets

Arthroscopy Sets - Sterile, Disposable

Product	Product Code	Length	Tip	Spike	Box Qty
	DAS4031	2m	Silicone	Single	20
	DAS4032	2.4m	Silicone	Double	20
	DAS4034	2m	Luer Lock	Single	20
	DAS4033	2.4m	Luer Lock	Double	20
	DAS4202	3.7m	Silicone	Double	15



Corporate Office
Fairmont Medical
11 Scoresby Road Bayswater 3153 VIC
Ph 1300 972 088
Ph International +613 9720 8840
Fx International +613 9237 9399

UK Head Office
Fairmont Medical Products Ltd
Elizabeth House
28 Baddow Road Chelmsford
Essex CM2 0DG
Ph: +44 (0)1245 206800

