


dif

CO₂ Insufflator Filters & Tubing Sets


Fairmont Medical offers a comprehensive range of CO₂ Insufflator Filter Sets to suit various insufflator models.



Fairmont Medical, for over 25 years, has been at the forefront of the design and manufacture of CO₂ filters, ensuring patient safety during minimal invasive surgery. With the mandated use of a CO₂ filter, the task of patient protection against infectious blood-borne pathogens prevents cross contamination between patients and provides protection for costly insufflator units.

- Universal compatibility with all popular insufflator units
- Easy to connect rotational luer lock system
- Disposable drapeable tubing set
- Bidirectional standard PTFE filter media offering both machine and patient protection
- Safe hydrophobic filter protection against the dangers of reflux fluids from patient into the insufflator unit
- Single use, sterile, ready to use
- Phthalate (DEHP) and latex free
- Quality Australian designed and manufactured

Standard CO₂ Filter



Hi-Flo CO₂ Filter











“Available with Hi-Flo CO₂ pleated filter.”







CO₂ Insufflator Filters & Tubing Sets

CO ₂ Insufflator Filters and Tubing Sets - Sterile, Disposable					
Product	Product Code	Insufflator Outlet Connector	Suits	Tubing Length	Box Qty
	DIF7002	40mm Silicone Universal Adaptor	6mm Barb	3m Tubing with rotating Male Luer Lock	50
	DIF7300	40mm Silicone Universal Adaptor	6mm Barb	3m Soft Tubing with rotating Male Luer Lock	50
	DIF7004	22mm Female Connector	22mm Male Outlet	3m Tubing with rotating Male Luer Lock	50
	DIF7005	Male Luer Lock Connector	Female Luer Outlet	3m Tubing with rotating Male Luer Lock	50
	DIF7007	Hexagonal Connector	Sopro Outlet	Y Extension and 3m Tubing with rotating Male Luer Lock	50
	DIF7009	Female Luer Lock	Male Luer Lock	3m with rotating Male Luer Lock	45
	DIF7205	15mm Male Connector	15mm Female Outlets	3m with rotating Male Luer Lock	50
	DIF7201	Universal Trumpet Connector (Red)	Fits a variety of Male Outlets	3m with rotating Male Luer Lock	40
	DIF7303	Disposable Insufflator Filter - High Flow with Pleated Filter, 15/22mm connector,	22M/15F Outlet	3.2m Tubing with rotating Male Luer Lock	10

CO ₂ Insufflator Filters - Sterile, Disposable				
Product	Product Code	Insufflator Outlet Connector	Suits	Box Qty
	DIF7001	40mm Silicone Universal Connector	6mm Barb Outlet	100
	DIF7003	Female ISO 22mm Connector	22mm Male Outlet	100

CO ₂ Insufflator Adaptors - Non-Sterile, Disposable			
Product	Product Code	Insufflator Outlet Connector	Box Qty
	DIF7015	15mm Insufflator Adaptor	1
	DIF7016	22mm Insufflator Adaptor	1



Corporate Office
Fairmont Medical
11 Scoresby Road Bayswater 3153 VIC
Ph 1300 972 088
Ph International +613 9720 8840
Fx International +613 9237 9399

UK Head Office
Fairmont Medical Products Ltd
Elizabeth House
28 Baddow Road Chelmsford
Essex CM2 0DG
Ph: +44 (0)1245 206800

