

# Disposable Camera Covers

Fairmont Medical pioneered the enclosed Camera Cover which is still the only 100% leak free sterile barrier in the world.

***“The Ultimate Sterile Barrier.”***



For over 25 years Fairmont Medical has been at the cutting edge of camera cover technology. Our draping system offers 100% true sterility where others fail. Our unique adaptor concept with locking mechanism offers safe and quick telescope exchange capability, mid surgical procedure.

- Simple setup procedure, telescopic folds allow easy draping
- Cover material is easily shaped over camera leads
- Enclosed design provides 100% protection
- Locking system avoids scope rotation
- Saves valuable operating theatre time with quick telescope exchange
- Fast, between case, camera turn around times
- Reduces the need for expensive camera system maintenance and inventory
- Single use, sterile, ready to use
- Quality Australian designed and manufactured



***“Accepted & respected around the world.”***

## Camera Cover - Sterile, Disposable

Product	Product Code	Style	Suit	Size	Box Qty
	DCC8007	Enclosed Adaptor with Telescopic Lock (Telescopically Folded)	All Operating Camera Systems	2850mm x 140mm	40



# Disposable Camera Covers




For over 25 years Fairmont Medical has been at the cutting edge of camera cover technology. Our draping system offers 100% true sterility where others fail. Our unique adaptor concept with locking mechanism offers a safe and quick telescope exchange capability, mid surgical procedure.

Fairmont Medical offers a comprehensive range of camera covers to suit all needs.

***“Additional camera covers to suit all your needs.”***

***“Fairmont Medical offers a comprehensive, cost effective range of camera covers to meet all design requirements.”***

- Simple setup procedure, telescopic folds allows easy draping
- Cover material is easily shaped over camera leads
- Single use, sterile, ready to use
- Quality Australian designed and manufactured

Camera Cover - Sterile, Disposable					
Product	Product Code	Style	Suit	Size	Box Qty
	DCC8000	Perforated Tip (Telescopically Folded)	All Operating Camera systems	1850mm x 140mm	40
	DCC8005	Enclosed Tip (Telescopically Folded) with vinyl end and Lens	All Operating Camera Systems	2100mm x 140mm and vinyl sock 220mm	40
	DCC8006	Enclosed Tip (Telescopically Folded) with vinyl end - without Lens	All Operating Camera Systems	2100mm x 140mm and vinyl sock 220mm	40



Corporate Office  
Fairmont Medical  
11 Scoresby Road Bayswater 3153 VIC  
Ph 1300 972 088  
Ph International +613 9720 8840  
Fx International +613 9237 9399

UK Head Office  
Fairmont Medical Products Ltd  
Elizabeth House  
28 Baddow Road Chelmsford  
Essex CM2 0DG  
Ph: +44 (0)1245 206800

