


dusflo

TURP/Cystoscopy Sets

For decades now, Fairmont Medical has been a leading manufacturer and supplier of an extensive range of urology sets to both private and public hospitals throughout the world.

Dusflo was the first urology set introduced to the Australian market to be latex free. Our large lumen tubing and drip chamber ensures high flow rates of irrigation fluids greater than 800ml/min which assists in efficient, rapid flushing of tissue debris from the bladder during TURP procedures. A soft non latex connector provides easy resectoscope connection, which in turn, allows easy attachment to the foley catheter post operatively.

Fairmont Medical's range of urology sets are a cost effective choice for both specialised TURP and cystoscopy procedures.



- Extra length tubing allows for smooth patient transfer from operating table to patient bed
- Easy instrument rotation and fitment via flexible non latex connector
- Time saving rapid irrigation flushing capabilities
- Easy to view and prime malleable drip chamber
- Fluid volume control via roller clamp
- Available in both single and double bag spike configurations
- Single use, sterile, ready to use
- Phthalate (DEHP) and latex free
- Quality Australian designed and manufactured





TURP/Cystoscopy Sets

TURP/Cystoscopy - Sterile, Disposable					
Product	Product Code	Length	Connector	Bag Spike	Box Qty
	DUS4005	1.5m	TPR	Single	20
	DUS4015	2.5m	TPR	Single	20
	DUS4008	1.8m	TPR	Double	20
	DUS4018	2.8m	TPR	Double	20
	DUS4021	2.8m	Luer Lock	Double	20



Corporate Office
Fairmont Medical
11 Scoresby Road Bayswater 3153 VIC
Ph 1300 972 088
Ph International +613 9720 8840
Fx International +613 9237 9399

UK Head Office
Fairmont Medical Products Ltd
Elizabeth House
28 Baddow Road Chelmsford
Essex CM2 0DG
Ph: +44 (0)1245 206800

