

dusflo plus TURP/Cystoscopy Sets

Fairmont Medical's new and improved range of urological irrigation sets offer a improved cost effective solution.


As part of Fairmont Medical's continued commitment to manufacture and design innovative medical devices, we are proud to introduce **Dusflo plus**.

Dusflo plus offers a new range of advanced TURP and cystoscopy sets, which further enhance urological procedures.

Quickflow

Quickflow combines all of the features common to Dusflo standard with following additional benefits.


- Large volume drip chamber incorporating bubble trap system
- Available in both single and double bag spike configurations
- Single use, sterile, ready to use
- Phthalate (DEHP) and latex free
- Quality Australian designed and manufactured



Eziflow

Eziflow combines all of the features common to Dusflo standard with following additional benefits.

- Large volume drip chamber with or without bubble trap system
- Shielded bag spike cover provides easy aseptic attachment
- Available in both single and double bag spike configurations
- Single use, sterile, ready to use
- Phthalate (DEHP) and latex free
- Quality Australian designed and manufactured





***dusflo** plus "Saves valuable surgical operating room time."*





TURP/Cystoscopy Sets



Eziflow TURP/Cystoscopy Sets - Sterile, Disposable

Product	Product Code	Length	Connector	Bag Spike	Box Qty
	DUS4205	1.5m	TPR	Single	20
	DUS4215	2.5m	TPR	Single	20
	DUS4208	1.8m	TPR	Double	20
	DUS4218	2.8m	TPR	Double	20

Eziflow with Bubble Trap TURP/Cystoscopy Sets - Sterile, Disposable

Product	Product Code	Length	Connector	Bag Spike	Box Qty
	DUS4505	1.5m	TPR	Single	20
	DUS4515	2.5m	TPR	Single	20
	DUS4508	1.8m	TPR	Double	20
	DUS4518	2.8m	TPR	Double	20

Quickflow TURP/Cystoscopy Sets - Sterile, Disposable

Product	Product Code	Length	Connector	Bag Spike	Box Qty
	DUS4305	1.5m	TPR	Single	20
	DUS4315	2.5m	TPR	Single	20
	DUS4308	1.8m	TPR	Double	20
	DUS4318	2.8m	TPR	Double	20


Fairmont
 MEDICAL

Corporate Office
 Fairmont Medical
 11 Scoresby Road Bayswater 3153 VIC
 Ph 1300 972 088
 Ph International +613 9720 8840
 Fx International +613 9237 9399

UK Head Office
 Fairmont Medical Products Ltd
 Elizabeth House
 28 Baddow Road Chelmsford
 Essex CM2 0DG
 Ph: +44 (0)1245 206800

