

**endoselect**


# Disposable Laparoscopic Scissors

The use of reusable laparoscopic scissors requires costly reprocessing and sterilisation, as well as frequent sharpening.


Fairmont Medical's, **Endoselect** Disposable Metzenbaum Scissors are designed to allow surgeons maximum sharpness during laparoscopic procedures when it counts. Supplied sterile ready to use, our range of laparoscopic scissors offers choice of jaw and shaft lengths to meet all laparoscopic needs.

- Balanced ergonomic handle designed for surgical control
- Safe state-of-the-art insulated shaft and jaw housing tip provides the ultimate diathermy current containment
- Curved blades allow for cutting delicate tissue whilst offering the capability of blunt dissection and coagulation
- Easy to identify colour coding
- Single use, sterile, ready to use
- Australian designed quality

Fairmont Medical's disposable Laparoscopic Scissors are designed for the precise control and handling during minimal invasive surgery.



**endoselect** "Guaranteed sharpness when it counts."

Sterile, Disposable					
Product	Product Code	Colour	Description	Blade Type/Size	Box Qty
	DLI4201	BLUE	Metzenbaum Laparoscopic Scissor 5mm x 330mm	Curved / 18mm	5
	DLI4202	GREEN	Metzenbaum Laparoscopic Scissor 5mm x 330mm	Curved / 12mm	5
	DLI4203	PURPLE	Metzenbaum Laparoscopic Scissor 5mm x 450mm	Curved / 18mm	5

**Fairmont**  
MEDICAL

Corporate Office  
Fairmont Medical  
11 Scoresby Road Bayswater 3153 VIC  
Ph 1300 972 088  
Ph International +613 9720 8840  
Fx International +613 9237 9399

UK Head Office  
Fairmont Medical Products Ltd  
Elizabeth House  
28 Baddow Road Chelmsford  
Essex CM2 0DG  
Ph: +44 (0)1245 206800

