

**endoselect**

# Laparoscopic Diathermy Hook

Fairmont Medical offers a range of disposable sterile 5mm insulated right angle hooks, presenting surgeons with a safe, ready to use choice. Extensively used throughout minimal invasive surgery, our range of diathermy hooks are insulated utilising a high tech coating of PTFE, polyolefin insulation providing peace of mind, safe use.

Fairmont Medical's range of Laparoscopic Diathermy Hooks offer a practical and cost effective solution.

## Standard Right Angle Hook

- Comfortable ergonomic handle design
- Variety of tip selections
- Offered with or without cable attachment
- Drapeable cable design on offer
- Ready to use, sterile single use
- Australian designed quality



**Standard Hook  
8mm Exposed**



**Insulated Hook  
5mm Exposed**



**Spatula**

## Switch Control Right Angle Hook

- Finger switch control avoids the use of awkward foot pedals
- Comfortable ergonomic handle design
- Easy to identify cut and coagulate colour coded buttons
- Safe reinforced 3 pin plug
- Convenient drapeable 3 metre cable
- Single use, sterile, ready to use
- Australian designed quality



**New Finger  
Switch Control**

**endoselect** "Guaranteed performance when it counts."



**30** YEARS OF  
AUSTRALIAN  
MANUFACTURING

**EXPORTED  
WORLDWIDE**



**Fairmont**  
MEDICAL

# Laparoscopic Diathermy Hook


## Standard Right Angle Hook - Sterile, Disposable

Product	Product Code	Connector	Probe Length	Tip Length	Insulation Type	Box Qty
	DHE4003	Banana Plug	330mm	8mm	5mm Extended	10
	DHE4005	Banana Plug	330mm	8mm	8mm Standard	10
	DHE4023	Split Stem Plug	330mm	8mm	5mm Extended	10
	DHE4025	Split Stem Plug	330mm	8mm	8mm Standard	10


## Standard Right Angle Hook With Integrated 3m Insulated Cable - Sterile, Disposable

Product	Product Code	Connector	Probe Length	Tip Length	Insulation Type	Box Qty
	DHE4103	8mm Generator Plug	330mm	8mm	5mm Extended	10
	DHE4105	8mm Generator Plug	330mm	8mm	8mm Standard	10


## Standard Right Angle Hook With Separate 3m Insulated Cable - Sterile, Disposable

Product	Product Code	Connector	Probe Length	Tip Length	Insulation Type	Box Qty
	DHE4113	8mm Generator Plug	330mm	8mm	5mm Extended	10


## High Frequency 3m Insulated Cable - Sterile, Disposable

Product	Product Code	Connector	Box Qty
	DHE4100	8mm Generator Plug	50

## Electrode Spatula with Banana Plug - Sterile, Disposable

Product	Product Code	Connector	Probe Length	Tip	Box Qty
	DHE4401	Banana Plug	330mm	Spatula	10

## Switch Control Right Angle Hook - Sterile, Disposable

Product	Product Code	Description	Probe Length	Tip	Box Qty
	DHE4505	Finger switch diathermy pen with removable right angle hook electrode. Connector - 3 metre cord with 3 prong plug	300mm	5mm Extended	10