



# Disposable Curettage Tubing Set

Fairmont Medical offers a cost effective range of Curettage Tubing Sets and Karman curettes.

Fairmont Medical, for over 25 years, has been at the forefront of offering innovative surgical devices. Our easy to use curettage tubing sets, including our extensive range of Karman curettes, offers a cost effective solution.

**Rovac** curettage tubing sets feature a ergonomically designed rotating hand piece, offering functionality and precision. An inbuilt on/off sliding mechanism allows control of suction flow rate.

- Large bore 3 metre tubing allows trouble free flow rate
- Clear view rotatable handle with on/off slider
- Flexible easy to fit suction cannister connector
- Comprehensive range of colour coded Karman curettes
- Easy to identify colour coded adaptors
- An extensive full range of curettes to choose from
- Phthalate (DEHP) and latex free
- Single use, sterile, ready to use
- Quality Australian designed and manufactured






**rovac** *"Saves valuable surgical operating room time."*



**Fairmont**  
MEDICAL



# Disposable Curettage Tubing Set

Curettage Suction Tubing - Sterile, Disposable					
Product	Product Code	Length	Tube ID	Tube OD	Box Qty
	DST3100 (Rovac Handle and Connector)	3m	9.8mm	15mm	10
	DST3101 (Rovac Handle)	3m	9.8mm	15mm	10
	DST3530 (Two Connectors)	3m	9.8mm	15mm	10

Karman Curettes - Sterile, Disposable				
Product	Product Code	Colour Tip	Diameter	Box Qty
	DKC3006	Grey Tip	6mm	10
	DKC3007	Yellow Tip	7mm	10
	DKC3008	Orange Tip	8mm	10
	DKC3009	Red Tip	9mm	10
	DKC3010	Green Tip	10mm	10
	DKC3011	Cyan Tip	11mm	10
	DKC3012	Blue Tip	12mm	10



Corporate Office  
Fairmont Medical  
11 Scoresby Road Bayswater 3153 VIC  
Ph 1300 972 088  
Ph International +613 9720 8840  
Fx International +613 9237 9399

UK Head Office  
Fairmont Medical Products Ltd  
Elizabeth House  
28 Baddow Road Chelmsford  
Essex CM2 0DG  
Ph: +44 (0)1245 206800

