

STIHLER ELECTRONIC

Astoflo Plus Eco - Blood & Infusion Warmer

Astoflo Plus eco provides the easiest and most economical way to effectively warm all infusion and transfusions.



Astoflo Plus Eco is the universal warming device for all infusions and blood products for operating theatres, intensive care units, and all other medical departments. All medical fluids are warmed in the electric heating profile - without special disposables being required. Each unit requires a heating profile which are available in different lengths and diametres. All medical fluids are warmed in the electric heating profile - without special disposables being required - in the standard infusion set and kept warm right up to the patient cannula.



Astoflo Plus Eco - Blood & Infusion Warmer

Product Code	Temperature Range	Initial Warm Up Time	Length of Connection Cable	Required Heating Profile Accessory
SEAFP300AU	33-43 degrees	Approx 3 minutes	40cm	SEWP31 - 180cm length* SEWP32 - 240cm length* SEWP33 - 180cm length** SEWP34 - 240cm length**
SEAFP302AU	33-43 degrees	Approx 3 minutes	80cm	

* For standard infusion sets with external diametres 4-5mm.

** For infusion sets with external diametres 6-7mm.

Fairmont
MEDICAL

Corporate Office
Fairmont Medical
11 Scoresby Road Bayswater 3153 VIC
Ph 1300 972 088
Ph International +613 9720 8840
Fx International +613 9237 9399

UK Head Office
Fairmont Medical Products Ltd
Elizabeth House
28 Baddow Road Chelmsford
Essex CM2 0DG
Ph: +44 (0)1245 206800

