

Astopad System

Astopad System improves safety against negative side effects of operations and contributes to the general well-being of the patient.



Heating Blanket

Under patient or over patient blanket




Heating Mattress

Waterproof and X-ray permeable



Astopad System can be used in pre-, peri- and post-operative areas and serves as a prophylaxis against patient cooling and as a hypothermia therapy. The Astopad system can control two Astopad blankets or mattresses independently of each other.

Astopad System							
Product	Product Code	Size & Weight	Temperature Selector	Operating Mode	Output Voltage	Alarm Functions	Two Heating Elements
	SEDU0310	250mm x 200mm x 185mm 4.2kg	32°C - 39°C at 0.5°C increments	Continuous (24h/day)	24V	Excess/Insufficient Temperature Cable Fracture	Controlled Independently

Astopad System

Astopad Heating Blankets

Astopad Heating Blankets safely and quickly warms the patient during surgery. The heat is generated by X-Ray permeable carbon fibre braid. Eight temperature sensors are arranged in the heated blanket close to the patient. These guarantee precise temperature control and monitoring in conjunction with alarm functions of the Astopad Control System. All of the Astopad Heated Covers can be washed down and disinfected. They can either be used as they are, or with washable reusable covers or with disposable covers.



Astopad Heating Blankets				
Ref	Product Code	Size	Accessory - Reusable Cover	Accessory - Disposable Cover
A	SECOV070 Paediatric	680mm x 480mm	SECOV40070	NA
B	SECOV105 Urology/Gynaecology	1050mm x 500mm	SECOV40105	NA
B	SECOV150 Full Size	1500m x 500mm	SECOV40150	SECOV30150 (BOX OF 50)
C	SECOV155 Arm/Chest	1500mm x 500mm	SECOV40155	SECOV30155 (BOX OF 50)
D	SECOV180 Large Full Body	1800mm x 800mm	SECOV40180	SECOV30180 (BOX OF 50)
E	SECOV235 Extra Full Body	2050mm x 1350mm	SECOV40235	NA



Astopad Heating Mattresses

Astopad Heating Mattresses provide heat and a unique combination of prophylaxis against hypothermia and intra-operative decubitus. The core of Astopad Heating Mattresses is a visco-elastic foam combined with a dynamic foam basis for best possible prophylaxis against decubitus on the OR table. The integrated heating element is regulated by the control unit and warms the patient from below. The Astopad Heating Mattresses are customised exactly to the customers needs. For the protection of the Astopad Mattresses, a cover is offered which is antistatic, disinfected and can be washed at 95°C.

Astopad Heating Mattresses			
Product Code	Custom Length Between	Custom Width Between	Custom Height Between
SESOF7	600mm to 800mm	450mm to 600mm	40mm to 100mm
SESOF4	810mm to 1300mm	480mm to 600mm	40mm to 100mm
SESOF5	1310mm to 1700mm	480mm to 600mm	40mm to 100mm
SESOF2	1710mm to 2300mm	480mm to 600mm	40mm to 100mm