

STIHLER ELECTRONIC

Astotherm Plus - Blood & Infusion Warmer

Astotherm Plus high-performance blood and infusion warmers covers the full range of clinically relevant flow rate.

Astotherm Plus warms all infusions and transfusions by the dry flow heating method efficiently, safely and economically in daily use. Astotherm Plus is ready for use within 2 minutes. Warming close to the patient means that infusions and blood products no longer require warming in advance. Blood products can be kept in the refrigerator until they are needed.



Fairmont
MEDICAL

Corporate Office
Fairmont Medical
11 Scoresby Road Bayswater 3153 VIC
Ph 1300 972 088
Ph International +613 9720 8840
Fx International +613 9237 9399

UK Head Office
Fairmont Medical Products Ltd
Elizabeth House
28 Baddow Road Chelmsford
Essex CM2 0DG
Ph: +44 (0)1245 206800





Astotherm Plus - Blood & Infusion Warmer

Astotherm Plus 220/220S

The universal device - effective performance for rapid routine use.

Astotherm Plus 220 is ready for use in under 60 seconds, with the result that every patient in theatre or in intensive care can be rapidly supplied with warm transfusions and infusions.



Astotherm Plus 220 covers over 97% of clinically-relevant transfusion and infusion rates

Product	Product Code	Size & Weight	Temperature Selections	Preheat Time	Groove Size	Power Consumption	Operating Mode
	SEAP220AUS	145mm x 135mm x 295mm 2.9kg	39°C 41°C 43°C	< 2 Min.	4mm, 11 Windings	450W	Continuous (24h/day)
	SEAP220SAUS	145mm x 135mm x 295mm 2.9kg	39°C 41°C 43°C	< 2 Min.	4mm, 11 Windings *Comes with heat protective sleeve and Astoline active insulation	450W	Continuous (24h/day)



Astotherm Plus 260/260S

The special device for the fastest flow rate.

In contrast to the Astotherm Plus 220, the Astotherm Plus 260 takes thicker tubes with an external diameter of 6.8mm but it is identical in terms of technology and potential applications. It is designed for warming blood and other fluids in dialysis, haemofiltration and urology.

Product	Product Code	Size & Weight	Temperature Selections	Preheat Time	Groove Size	Power Consumption	Operating Mode
	SEAP260AUS	145mm x 135mm x 295mm 2.9kg	39°C 41°C 43°C	< 2 Min.	6.8mm 11 Windings	450W	Continuous (24h/day)
	SEAP260SAUS	145mm x 135mm x 295mm 2.9kg	39°C 41°C 43°C	< 2 Min.	6.8mm 11 Windings *Comes with heat protective sleeve and Astoline active insulation	450W	Continuous (24h/day)

AstoTube - PVC Extension Line

Product	Product Code	Suits	Box Qty
	SEIFT30460	SEAP220AUS / 220SAUS	50
	SEIFT30410	SEAP260AUS / 260SAUS	25