

ultracover

Microscope Covers

Fairmont Medical's Microscope Covers are available in multiple lengths, sizes and variations.

“Preventing patient to patient cross contamination”



Fairmont Medical's Microscope Covers include convenient adjustable ties. This enables operating microscopes to be securely, neatly covered aseptically. The unique one piece adhesive tab system allows for easy adjustment and repositioning of the cover, if and when required.

- Fast secure ring fitment to microscope optics
- Low noise rustle free material cover
- Particle and angel hair free thermosealed shaped cover
- Clean cut ocular perforation
- Adjustable ocular secura-tabs
- Full range of drapes to suit most popular brands of microscopes
- Single use, sterile, ready to use
- Quality Australian designed and manufactured



ultracover *“The ultimate sterile barrier.”*



Fairmont
MEDICAL

Microscope Covers

Microscope Cover 1 Viewing Point - Sterile, Disposable

Product	Product Code	Length	Width	Lens Size	Lens Ring Colour	# Securing Ribbons	Box Qty
	DMC8011	1800mm	700mm	47mm	Orange	2	10
	DMC8014	1800mm	700mm	65mm	Blue	2	10
	DMC8111	1800mm	700mm	49mm	Yellow	2	10

Microscope Cover 2 Viewing Points - Sterile, Disposable

Product	Product Code	Length	Width	Lens Size	Lens Ring Colour	# Securing Ribbons	Box Qty
	DMC8051	1500mm	1040mm	47mm	Orange	2	10
	DMC8151	1500mm	1040mm	49mm	Yellow	2	10

Microscope Cover 3 Viewing Points - Sterile, Disposable

Product	Product Code	Length	Width	Lens Size	Lens Ring Colour	# Securing Ribbons	Box Qty
	DMC8010	1800mm	1040mm	47mm	Orange	2	10
	DMC8013	1500mm	1040mm	65mm	Blue	2	10
	DMC8020	1800mm	1040mm	69mm	Pink	2	10
	DMC8110	1800mm	1040mm	49mm	Yellow	2	10
	DMC8030	3000mm	1040mm	47mm	Orange	3	5
	DMC8032	3000mm	1150mm	65mm	Blue	3	5
	DMC8040	3000mm	1040mm	69mm	Pink	3	5
	DMC8130	3000mm	1040mm	49mm	Yellow	2	5
	DMC8060	3500mm	1040mm	47mm	Orange	4	5
	DMC8081	4500mm	1150mm	69mm	Pink	5	5
	DMC8160	3500mm	1040mm	49mm	Yellow	4	5

Microscope Cover - Shower Cap Style - Sterile, Disposable

Product	Product Code	Length	Width	Lens Size	Lens Ring Colour	Box Qty
	DMC8066	700mm	450mm	65mm	Blue	20
	DMC8070	700mm	450mm	47mm	Orange	20
	DMC8170	700mm	450mm	49mm	Yellow	20
	DMC8080	700mm	450mm	69mm	Pink	20

Speciality Microscope Cover - BIOM - Sterile, Disposable

	DMC8067 (with securing screw)	700mm	450mm	65mm	Blue	20
--	----------------------------------	-------	-------	------	------	----

Mini Style with No Lens - Sterile, Disposable

	DMC8065	700mm	450mm	65mm	Blue	20
--	---------	-------	-------	------	------	----