

ultracover

Ultrasound Probe Covers

Fairmont Medical's Ultrasound Probe Covers provide the ultimate sterile barrier preventing patient to patient cross contamination.

“Preventing patient to patient cross contamination.”



Fairmont Medical offers a range of Ultrasound Probe Cover sizes to suit most applications. With recent directives regarding the use of ultrasound probes, the necessity of a disposable ultrasound probe cover is becoming more of a requirement during intra-operative surgical procedures. Manufactured from a durable non latex material, ensures a sterile barrier without compromise.

- Soft conforming PU material utilised
- Variety of size and lengths to suit all probe designs
- Sterile ultrasound gel sachet (DUC8133)
- Non latex rubberbands included
- Single use, sterile, ready to use
- Quality Australian designed and manufactured



ultracover *“The ultimate sterile barrier.”*



-Fairmont
MEDICAL

Ultrasound Probe Covers

Ultrasound Probe Covers - Sterile, Disposable					
Product	Product Code	Type	Length	Width	Box Qty
	DUC8016	PP Telescopically Folded	1700mm	140mm	40
	DUC8132	Non-Latex Elastic bands - TPU Tapered	640mm	45mm/100mm	40
	DUC8134	Non-Latex Elastic bands - TPU Tapered	250mm	40mm/80mm	80
	DUC8135	Non-Latex Elastic bands - TPU Tapered	220mm	100mm/140mm	40
	DUC8137	Non-Latex Elastic bands - TPU Tapered	400mm	45mm/100mm	40
	DUC8133	Non-Latex Elastic bands - TPU Telescopically Folded with Gel	1230mm	130mm (Square)	40
	DUC8136	PP/TPU Tapered - Telescopically Folded	1100mm	45mm/100mm	40
	DUC8138	PP/TPU Tapered - Telescopically Folded	2340mm	100mm/140mm	40
	DUC8139	PP/TPU Tapered - Telescopically Folded	1500mm	100mm/140mm	40