

Disposable Speciality Suction

Fairmont Medical has a comprehensive range of high quality specialty suction tips to suit all surgical applications. Poole, Frazier and Paediatric Bubble Sucker.

Universal Poole Sucker

Uniquely designed to evacuate pooled blood, fluids and tissue debris during abdominal surgeries and cesarean section. Our single use Poole Sucker suction device offers a cost effective, efficient solution.

Paediatric Bubble Suckers

Fairmont Medical's clear Paediatric Bubble Suckers enable easy visualisation of suctioned fluids. Atraumatic tip design allows for the suction of delicate tissue, incorporating in a smooth finish for comfortable easy handling.

Frazier Suckers

Fairmont's Frazier suckers are used in surgical procedures where delicate fine suctioning is required to evacuate tissue debris from the site. The included obturator facilitates the shaping and clearing of the flake resistant malleable metal shaft.

Fairmont Medical has a comprehensive range of Speciality Suction Tips to suit all applications.



"Cost effective solutions."


Poole Sucker - Sterile, Disposable

Product	Product Code	Description	Box Qty
	DYS6830	Universal Poole Sucker	50

Frazier Sucker - Sterile, Disposable

Product	Product Code	Description	Box Qty
	DFS0008	Long Handle 8 Fr	25
	DFS0010	Long Handle 10 Fr	25
	DFS0012	Long Handle 12 Fr	25
	DFS1006	Short Handle 6 Fr	25
	DFS1008	Short Handle 8 Fr	25
	DFS1010	Short Handle 10 Fr	25
	DFS1012	Short Handle 12 Fr	25

Speciality Yankauer Suckers - Sterile, Disposable

Product	Product Code	Description	Box Qty
	DYS6318	Paediatric Bubble Yankauer Sucker Non Vented 18 FR gauge	100