



# Disposable Yankauer Suckers

Fairmont Medical offers an extensive range of suction devices to suit all surgical and general needs.

Fairmont's Australian made range of yankauer suction handles are designed for comfort and reliability. **Vacusite's** durability and strength make it an ideal choice for suctioning in a range of surgical procedures, accident and emergency, ICU, anaesthesia and ward use.

- Ergonomically designed hand control, ensures accurate, comfortable use
- Ridged barbed connector ensures secure suction tubing fit
- Smooth tip ensures atraumatic suctioning
- Available in rose tip or anaesthetic tip
- Choice of vented or non vented version
- Securely packaged with double wrap if required
- On/off switch version on offer
- Single use, sterile, ready to use
- Quality Australian designed and manufactured



**Straight Tip**



**Rose Tip**



## On/Off Yankauer

Our uniquely designed on/off yankauer sucker reduces suction noise when not in use. Suction vacuum is easily controlled allowing delicate utilisation.

Yankauer Suckers - Sterile, Disposable					
Product	Product Code	Tip	Vent	Wrap	Box Qty
	DYS6822	Straight	Non Vented	Single	100
	DYS6828	Straight	Non Vented	Double	100
	DYS6823	Straight	Vented	Double	100
	DYS6829	Straight	Vented	Single	100
	DYS6824	Rose	Non Vented	Single	100
	DYS6825	Rose	Non Vented	Double	100
	DYS6826	Rose	Vented	Double	100
	DYS6827	Rose	Vented	Single	100
Speciality Yankauer Suckers - Sterile, Disposable					
Product	Product Code	Description			Box Qty
	DYS6840	On/ Off - Straight Tip Double Wrapped			30
	DYS6850	On/ Off - Rose Tip Double Wrapped			30



Corporate Office  
Fairmont Medical  
11 Scoresby Road Bayswater 3153 VIC  
Ph 1300 972 088  
Ph International +613 9720 8840  
Fx International +613 9237 9399

UK Head Office  
Fairmont Medical Products Ltd  
Elizabeth House  
28 Baddow Road Chelmsford  
Essex CM2 0DG  
Ph: +44 (0)1245 206800

