

Diathermy Switch Pen

Fairmont Medical's Finger Switch Diathermy Pen is designed for single use monopolar electrosurgery.

Standard

Fairmont Medical provides a diathermy switch pen, and comprehensive range of accessories, that meet with the highest quality and safety demands of today. Our standard finger switch pen comes complete with a selection of electrode types, ranging from insulated, teflon coated through to the standard stainless steel electrode.

- Easy lightweight slimline design allows exceptional handling and control
- Drapeable cable with reinforced 3 pin plug
- Easy to identify colour coded cut and coag buttons
- Simple interchangeable electrode fitment
- Extensive range of electrode types and lengths on offer: bladed, needle, ball, D shape and LETZ loop electrodes
- Single use, sterile, ready to use
- Australian designed quality



Various non-stick Teflon coated electrodes available



Blade electrode



Needle electrode tip




Ball electrode



Loop electrode tip

Diathermy Switch Pen - Sterile, Disposable

Product	Product Code	Description	Box Qty
	DSP1110	Disposable Finger Switch Diathermy Pen with 70mm Stainless Steel Blade electrode - 3mt cable & 3 prong plug	50
	DSP1140	Disposable Finger Switch Diathermy Pen with 70mm Teflon Blade electrode - Extended Insulation, 3mt cable & 3 prong plug	50

Fairmont
MEDICAL

Corporate Office
Fairmont Medical
11 Scoresby Road Bayswater 3153 VIC
Ph 1300 972 088
Ph International +613 9720 8840
Fx International +613 9237 9399

UK Head Office
Fairmont Medical Products Ltd
Elizabeth House
28 Baddow Road Chelmsford
Essex CM2 0DG
Ph: +44 (0)1245 206800

