

endosac

Retrieval Pouches

Fairmont Medical's disposable Retrieval Devices provide an economical retrieval system for laparoscopy procedures.

Fairmont Medical's, **EndoSac** Retrieval Pouches are designed to provide cost effective measures that avoids surgical site infection complications post surgery, such as laparoscopic cholecystectomy. Minimises port site infection and carcinoma contamination seepage by isolating and containing specimens during extraction.

- Available in strong polyurethane, or Ripstop nylon, tear resistant, non-porous pouch material.
- Easy view transparent pouch for tissue visualisation
- Ensuring quick, simple insertion via 10mm introducer
- Single use, sterile, ready to use
- Australian designed quality



Tri Handle

Uni Wand

Ripstop Pouch



PU Basic Pouch

endosac "Secure and easy deployment retrieval pouches."

30 YEARS OF AUSTRALIAN MANUFACTURING

EXPORTED WORLDWIDE

Fairmont MEDICAL

Retrieval Pouches

EndoSac Tri handle - Sterile, Disposable			
Product	Product Code	Description	Box Qty
	DRB0100	100mL bag (60mm x 150mm) 10mm x 330mm Butterfly handle introducer	10
	DRB0200	200mL bag (60mm x 180mm) 10mm x 330mm Butterfly handle introducer	10
	DRB0400	400mL bag (100mm x 140mm) 10mm x 330mm Butterfly handle introducer	10
	DRB0700	700mL bag (100mm x 200mm) 10mm x 330mm Butterfly handle introducer	10
	DRB1200	1200mL bag (100mm x 240mm) 10mm x 330mm Butterfly handle introducer	10
	DRB1600	1600mL bag (100mm x 290mm) 12mm x 330mm Butterfly handle introducer	10
	DRB3100	100mL bag Ripstop Pouch 10mm x 330mm Butterfly handle introducer	10
	DRB3200	200mL bag Ripstop Pouch 10mm x 330mm Butterfly handle introducer	10

Endosac UniWand - Sterile, Disposable			
	DEB0200	200mL bag with memory wire (100mm x 130mm) 10mm introducer	10
	DEB0400	400mL bag with memory wire (160mm x 140mm) 10mm introducer	10
	DEB0700	700mL bag with memory wire (170mm x 200mm) 10mm introducer	10
	DEB1600	1600mL bag with memory wire (170mm x 290mm) 12mm introducer	10
	DEB3000	3000mL bag with memory wire (210mm x 300mm) 12mm introducer	10

EndoSac PU Basic Pouch - Sterile, Disposable			
Product	Product Code	Description	Box Qty
	DEP0200	200mL, 10 x 13cm endo pouch with PE string	50
	DEP0400	400mL, 16 x 14cm endo pouch with PE string	50
	DEP0700	700mL, 17 x 20cm endo pouch with PE string	50
	DEP1600	1600mL, 17 x 29cm endo pouch with PE string	50
	DEP3000	3000mL, 21 x 30cm endo pouch with PE string	50

EndoSac Ripstop Pouch - Sterile, Disposable			
Product	Product Code	Description	Box Qty
	DRP0200	200mL, 10 x 13cm, Ripstop Pouch	50
	DRP0400	400mL, 16 x 14cm, Ripstop Pouch	50
	DRP0700	700mL, 17 x 20cm, Ripstop Pouch	50