

Airway Management Accessories

Fairmont Medical's Airway Management Accessory range provides some unique alternative products.

C-Mount Kits

C-Mount Kits include filter, collapsible tubing mount and detachable elbow



Standard Elbow *Luer Port Elbow*



Guedel Airways

- The smooth, soft finish of the airway reduces the risk of damage to oropharyngeal tissue
- Colour coded for easy visualisation of size choice
- Prevents potential dental damage and occlusion of airway in event of involuntary biting








Bite-Me Perioperative Bite Block


BiteMe™ is a purpose designed, air-filled, soft plastic bite block to prevent dental trauma and airway device (ETT or LMA) damage in anaesthetised patients in the operating theatre or post anaesthetic care unit.



- Soft malleable FDA approved plastic material
- Easy placement between teeth enabled by tapered distal end
- Allows conformity fitment to airway device via concave shape
- Easy to read depth markers allows accurate positioning
- Safe blunt distal tip decreases risk of oral trauma
- Easy insertion and removal due smooth surface
- Single use
- Highly visible whilst in place
- Patented design







Airway Management Accessories


Anaesthetic Filter Kits - Disposable			
Product	Product Code	Description	Box Qty
	DAF9027 (Non-Sterile)	Kompact Deltastatic Filter, collapsible C Mount 120mm-210mm 15F-22F and detachable elbow connector 15M-22M/15F NO luer port	50
	DAF9028 (Sterile)	Kompact Deltastatic HME Filter and detachable elbow connector 15M-22M/15F NO luer port	50
	DAF9018 (Sterile)	Deltastatic Filter, collapsible C Mount 120mm-210mm 15F-22F and detachable elbow 15M-22M/15F NO luer port	50
	DAF9019 (Sterile)	Deltastatic HME Filter, collapsible C Mount 120mm-210mm 15F-22F and detachable elbow 15M-22M/15F NO luer port	50
	DAF9045 (Sterile)	Isopleat Filter, collapsible C Mount 120mm-210mm 15F-22F and detachable elbow 15M-22M/15F NO luer port	50

Catheter Mount Kit - Sterile, Disposable			
Product	Product Code	Description	Box Qty
	DAF9290	120mm-210mm; 15F-22F connectors and detached elbow connector without luer port 22M/15F	50

Elbow Connector - Non-Sterile, Disposable			
Product	Product Code	Description	Box Qty
	DAF9091	Elbow connector with luer port 15M-22M/15F	50
	DAF9090	Elbow connector without luer port 15M-22M/15F	50

Catheter Mounts - Sterile, Disposable			
Product	Product Code	Description	Box Qty
	DAF9036	Collapsible Catheter Mount: 120mm-210mm 15M with attached elbow connector 22M/15F. No Luer port	50
	DAF9037	Collapsible Catheter Mount: 120 mm-210mm 15F-22F connector (Elbows not included)	50
	DAF9038	Collapsible Catheter Mount: 120mm-210mm 15M with attached elbow connector 22M/15F with Luer port	50
	DAF9039	Collapsible Catheter Mount: 120mm-210mm 22F with attached Elbow connector 22M/15F - No Luer Port	50

Guedel Airways - Non-Sterile, Disposable				
Product	Product Code	Colour	Tubing	Box Qty
	DGA9204NS	Pink	Guedel Airway Size 000 - 40mm	50
	DGA9205NS	Blue	Guedel Airway - Size 00 - 50mm	50
	DGA9206NS	Black	Guedel Airway - Size 0 - 60mm	50
	DGA9207NS	White	Guedel Airway - Size 1 - 70mm	50
	DGA9208NS	Green	Guedel Airway Size 2 - 80mm	50
	DGA9209NS	Yellow	Guedel Airway Size 3 - 90mm	50
	DGA9210NS	Red	Guedel Airway Size 4 - 100mm	50
	DGA9211NS	Light Blue	Guedel Airway Size 5 - 110mm	50
	DGA9212NS	Orange	Guedel Airway Size 6 - 120mm	50

BiteMe Perioperative Bite Block Sterile, Disposable			
Product	Product Code	Description	Box Qty
	BIT1000	Bite-Me Perioperative Bite Block	100



Corporate Office
Fairmont Medical
11 Scoresby Road Bayswater 3153 VIC
Ph 1300 972 088
Ph International +613 9720 8840
Fx International +613 9237 9399

UK Head Office
Fairmont Medical Products Ltd
Elizabeth House
28 Baddow Road Chelmsford
Essex CM2 0DG
Ph: +44 (0)1245 206800

