

endostream

Endoscopic Irrigation Set

Fairmont Medical's EndoStream system offers cost effective solutions.

Fairmont Medical's **Endostream** irrigation sets offers GI Endoscopy surgical teams a practical and cost effective disposable alternative to reusable systems.


Endostream is equipped with a back-flow check valve as a safeguard against potential infection control risks and work in conjunction with a variety of irrigation pumps.

- Reduces the risk of cross contamination
- Optimal flow and leakage protection from hermetically sealed cap
- Uniquely designed universal cap and connector fits all sterile branded water irrigation bottles and endoscopes
- Eliminates manual cleaning and reprocessing of reusable tubing
- Phthalate (DEHP) and latex free
- Single use, sterile, ready to use
- Quality Australian designed and manufactured



endostream "Saves valuable procedural time."

Sterile, Disposable

Product	Product Code	Description	Box Qty
	DES5013	3 metre Disposable Irrigation Circuit with Universal Cap	20

Fairmont
MEDICAL

Corporate Office
Fairmont Medical
11 Scoresby Road Bayswater 3153 VIC
Ph 1300 972 088
Ph International +613 9720 8840
Fx International +613 9237 9399

UK Head Office
Fairmont Medical Products Ltd
Elizabeth House
28 Baddow Road Chelmsford
Essex CM2 0DG
Ph: +44 (0)1245 206800

