

NONIN Reusable Flex Sensors



Nonin's Flex Sensor System combines the convenience and performance of a disposable sensor with the economics of a reusable sensor. The reusable Flex Sensor, made of a durable silicone, is secured with a matching disposable FlexiWrap adhesive. This hybrid sensor system is ideal for numerous applications including extended monitoring and sleep studies. Available in sizes to fit neonates to adults, the Flex Sensor System is a great value for patient monitoring. Each sensor is packaged with 25 FlexiWraps (additional packages of 25 may be purchased separately).

- Proven PureLight LEDs for accuracy at critical SpO₂ levels
- Durability for repeated use
- Optimal performance in motion situations
- Cost-effective convenience
- Simple cleaning for ease-of-use
- Trim shape for comfort and contoured fit

Available Models

Reusable Flex Sensors

NON-8000J	Flex Reusable Sensor - Adult, 1 metre cable with 25 FlexiWraps	NON-8001J	Flex Reusable Sensor - Neonate, 1 metre cable with 25 FlexiWraps
NON-8000J-3	Flex Reusable Sensor - Adult, 3 metre cable with 25 FlexiWraps	NON-8008J	Flex Reusable Sensor - Infant, 1 metre cable with 25 FlexiWraps
NON-8000J-WO2	Flex Reusable Sensor for Wrist Ox ₂ - Adult, 0.3 Metre cable & 25 FlexiWraps		

Disposable FlexiWrap

NON-8000JFW	FlexiWraps for Adult Flex Sensor (8000J) 25 Pack	NON-8008JFW	FlexiWraps for Infant Flex Sensor (8008J) 25 Pack
NON-8001JFW	FlexiWraps for Neonate Flex Sensor (8001J) 25 Pack		

Fairmont
MEDICAL

Corporate Office
Fairmont Medical
11 Scoresby Road Bayswater 3153 VIC
Ph 1300 972 088
Ph International +613 9720 8840
Fx International +613 9237 9399

UK Head Office
Fairmont Medical Products Ltd
Elizabeth House
28 Baddow Road Chelmsford
Essex CM2 0DG
Ph: +44 (0)1245 206800

