

Diathermy Smoke Pen

Fairmont Medical's Diathermy Smoke Pen optimises smoke evacuation contributing to a healthier O.R.




Protect your surgical team today and into the future with Fairmont's Pencavac Smoke Pen.

- Effortless, easy to extend hood and electrode reach
- Saving valuable operating time
- Super wide smoke evacuation coverage > 50mm
- Aerodynamic design affords high vacuum > 80L/min
- Easy clear view transparent graduated smoke evacuation nozzle
- Telescopic extension nozzle length up to 130mm
- Light, flexible tubing and swivel connector enables enhanced handling and control
- Selection of optional electrodes available
- Single use, sterile, ready to use
- Phthalate (DEHP) and latex free
- Australian designed quality



Stretchie Pen

- Easy to identify colour coded cut and coag buttons
- Compact streamlined fully integrated hose and cable

Pencavac Telescopic Diathermy Smoke Pen - Sterile, Disposable			
Product	Product Code	Description	Box Qty
	DSP3140	Disposable Smoke Evacuator Diathermy Pen with Extendable shaft - 63.5mm Teflon Blade electrode, 3.0m Tubing with 22mm Female/Male connector and 3m cable with 3-prong plug with holster	10
	DSP3142	Disposable Smoke Evacuator Diathermy Pen with Extendable shaft - 63.5mm Teflon Blade electrode, 3.0m Tubing with Suction connector and 3m cable with 3-prong plug with holster	10