

ABSORBEST



Fairmont
MEDICAL

ABSORBEST DryMax

DryMax super-absorbent pads improve hygiene and infection control within operating rooms.

No more slippery floors. No more soaked, cold and dripping cloths and sheets. May we introduce DryMax Operation, a range of super-absorbent mats and pads that efficiently manages excess fluids that occur on floors, beds and tables in operating theatres, emergency situations or during medical examinations.

DryMax super-absorbent products quickly absorb and encapsulate large amounts of water-based liquids, such as saline, urine and blood, including contaminated parts, without leaking. The result is a drier, safer and more hygienic working environment – with improved infection control – which also means less time and money spent on cleaning.

- Superior absorption – superior retention
- Comfortable, user-friendly – easy to apply
- No dripping – easy and safe to change and discard
- Fewer changes – less work and lower material costs
- Simplicity – easy to use



DryMax Sterile

Super-absorbent pads for use with the patient. Available in two sizes, with or without tape.






DryMax 2.4

Super-absorbent floor mat with high retention capacity, in a handy size.

DryMax XL

Anti-slip floor mat for collection of splashing liquids.

DryMax - Sterile, Disposable					
Product	Product Code	Water Absorption	Saline Absorption Rate	Dimensions	Box Qty
	AF50243 Maxi with Tape	6L	2.4L	720mmx 370mm	80
	AF50163 Midi with Tape	4L	1.6L	560mmx 370mm	100
DryMax - Non-Sterile, Disposable					
Product	Product Code	Water Absorption	Saline Absorption Rate	Dimensions	Box Qty
	AF22201 (Blue)	1.5 L	1.5L of Fluid	1000mm x 750mm	50
	AF20241 (Green)	9L	9L of Water 3.5L of Saline	720mm x 370mm	100



Corporate Office
Fairmont Medical
11 Scoresby Road Bayswater 3153 VIC
Ph 1300 972 088
Ph International +613 9720 8840
Fx International +613 9237 9399

UK Head Office
Fairmont Medical Products Ltd
Elizabeth House
28 Baddow Road Chelmsford
Essex CM2 0DG
Ph: +44 (0)1245 206800

