

# Ventilator Circuit Peep Valve


Fairmont Medical's Peep Valve offers a cost effective and safe product for critical care patients.

Fairmont Medical's Ventilator Peep Valve is a spring-loaded valve that a patient exhales against to allow respiration. Peep Valves are designed for use with manual resuscitators or ventilators. The Fairmont Medical Peep Valve is designed in a transparent material which allows monitoring of the patient's respiratory rate with measuring indicators for accurate flow rate. This Australian made device is a safe and disposable solution for emergency medical situations.

- Reduces the risk of cross contamination
- Optimal flow and leakage protection from hermetically sealed cap
- Transparent casing enables monitoring of patient's respiratory rate and blockage assessment
- Adjustable PEEP valve that connects to the one-way valve which incorporates a swivel connection to the patient to allow for maneuverability.
- Eliminates manual cleaning and reprocessing of reusable tubing
- Phthalate (DEHP) and latex free
- Single use, sterile, ready to use
- Quality Australian designed and manufactured



## Non-Sterile, Disposable

Product	Product Code	Description	Box Qty
	DPV9690	Disposable PEEP Valve	25

