



Ventilator Circuit

Fairmont Medical's Ventilator Circuit offers a cost effective and safe product for critical care patients.

Fairmont Medical's **Ventilator Circuit** offers corrugated flexible kink-resistant tubing which forms a closed system between the humidifier and ventilator. Fairmont's Ventilator Circuit connects the ventilator to the patient, as well as any devices that might be connected to the circuit. This disposable circuit offers medical staff a practical and cost effective disposable alternative to reusable systems.

- Reduces the risk of cross contamination
- Optimal flow and leakage protection from hermetically sealed cap
- Corrugated tubing with a pressure monitoring line that runs within the tube
- Adjustable PEEP valve that connects to the one-way valve which incorporates a swivel connection to the patient to allow for maneuverability.
- Eliminates manual cleaning and reprocessing of reusable tubing
- Phthalate (DEHP) and latex free
- Single use, sterile, ready to use
- Quality Australian designed and manufactured



| Non-Sterile, Disposable | | | |
|---|--------------|--|---------|
| Product | Product Code | Description | Box Qty |
|  | DVC9681 | Disposable Ventilator Circuit 1.5m - with Peep Valve | 10 |
|  | DVC9682 | Disposable Ventilator Circuit - Parallel with Peep Valve & Filters | 5 |

